

Qualitative Analysis of Salts

Subject: Chemistry | Grade: Grade 12

Objective:

To identify the cation and anion present in a given unknown salt sample through systematic qualitative analysis.

Materials:

Test Tubes and Rack
Bunsen Burner
Nichrome Wire
Dilute HCl, HNO₃
NaOH Solution
NH₄OH Solution
AgNO₃ Solution
BaCl₂ Solution
Limewater
Unknown Salt Samples

Procedure:

1. Preliminary Tests: Note color, solubility, and appearance.
2. Flame Test: Clean nichrome wire in HCl. Dip in salt paste. Hold in flame. Observe color.
(Yellow=Na, Lilac=K, Brick Red=Ca, Apple Green=Ba, Blue-Green=Cu)
3. Cation Analysis: Add NaOH dropwise then in excess. Repeat with NH₄OH. Observe precipitates.
4. Anion Analysis:
 - Carbonate test: Add dilute acid. Check for effervescence (CO₂) turning limewater milky.
 - Chloride test: Add dilute HNO₃ then AgNO₃. White ppt = Chloride.
 - Sulfate test: Add dilute HCl then BaCl₂. White ppt = Sulfate.

Discussion Questions:

1. Identify the cation if a white precipitate dissolves in excess NaOH.
2. Identify the anion that forms a yellow precipitate with AgNO₃.
3. Why is the nichrome wire cleaned with HCl?
4. Write the ionic equation for the test for Chloride ions.